


# HOW TO TEST YOUR DRUGS FOR *XYLAZINE*



Xylazine strips should not be used on cocaine.

If using alone, call  
**National Overdose Response Service**  
1-888-688-NORS (6677)

For more information, visit  
**timiskaminghu.com**



1. Put your drugs (5-10 mg or a tiny pinch) into 5 ml (teaspoon) of room temperature water.
2. Stir until the drug is completely dissolved.
3. Dip the strip into the cup for **15 seconds**, no deeper than the wavy line.
4. Wait 1 minute for **red lines** to appear. Do not read after 10 minutes.
5. Check the test results. **One line means that the substance has xylazine in it.**



**1 LINE:**  
**POSITIVE** for xylazine



**2 LINES:**  
**NEGATIVE** for xylazine

A negative test result does not mean the substance is safe. **Always use caution. Carry Naloxone.**